

Progress Towards A New Technique for Measuring Local Electric and Magnetic Field Fluctuations in High Temperature Plasmas

Galen Burke

R.J. Fonck, G. Mckee, G.R. Winz

APS DPP, Milwaukee

October 26, 2017




Introduction




Understanding Turbulence in Tokamaks is a Fundamental Challenge for Fusion Energy

- Plasma turbulence in tokamaks results in anomalous transport
 - Cross-field transport \gg neoclassical predictions
- Present plasma diagnostics measure key fluctuating parameters \tilde{n} , \tilde{T}_i , \tilde{T}_e , \tilde{v}
- Measurement of electrostatic field turbulence ($\tilde{E} \sim k_{\perp} \tilde{\phi}$) remains a challenge, gives \tilde{v}
 - $\tilde{E}_{\theta} \times B_{\phi} \cong \tilde{v}_r$: turbulent cross-field transport, $\tilde{E}_r \times B_{\phi} \cong \tilde{v}_{\theta}$: shear-flow and zonal flow dynamics
- Local magnetic field fluctuation (\tilde{B}) measurement also challenging, could provide critical information
 - Local \tilde{B} dynamics during edge harmonic oscillation (EHO)
 - 3D magnetic field perturbation penetration into plasma pedestal



Density perturbations and calculated velocimetry in DIII-D plasma




GYRO turbulence simulation

Motional Stark Effect Field Used as Carrier Signal for \tilde{E}

- Motional Stark Effect spectrum provides carrier line broadening for \tilde{E} :



$$\vec{E}_{total} = \vec{v}_b \times \vec{B} + \vec{E}_{plasma}$$

- Measure high-speed variations in π/σ line intensity ratio or in π -components linewidth to derive \tilde{E}
- Spatial Heterodyne Spectrometer (SHS) provides flexible analyzer of multiplet spectrum
- New CMOS imaging systems provide detection and DAQ



Multiple Techniques Used to Extract Components of $\vec{E}_{plasma}(t)$

- Midplane beam, sightline: linewidth changes
 - Radial beam: $\vec{E}_{MSE}\hat{z}$, \tilde{E}_z doesn't change α
 - $\tilde{\Delta\lambda}_{Stark} \propto \tilde{E}_z \rightarrow \tilde{v}_r$
- Midplane beam, off-midplane sightline: Intensity ratio change
 - $R = \frac{\sum I_\pi}{\sum I_\sigma} = \frac{\sin^2 \alpha}{1 + \cos^2 \alpha} (I_\pi / I_\sigma) (n_e) \equiv \frac{\sin^2 \alpha}{1 + \cos^2 \alpha} F$
 - $\tilde{R} = R \left[\frac{\partial \ln F}{\partial \ln n_e} \frac{\tilde{n}_e}{n_e} + \frac{4 \cos \alpha}{(1 + \cos^2 \alpha) \sin \alpha} \tilde{\alpha} \right]$
 - $\tilde{E}_r \sim \tilde{\alpha} E_{MSE} \rightarrow \tilde{v}_z$
- First emphasis on line width measurement: insensitive to density fluctuations, midplane view



Local Magnetic Field Fluctuations May be Measurable via Stark Multiplet

- Measurement of local magnetic field fluctuations (\tilde{B}) in high temperature plasmas is challenging
 - Provides information on: fast particle modes, island structures, plasma response to 3D RMP for ELM control, high- β turbulence
- Again use MSE field as carrier
- For midplane sightline and radial beam:
 - $\vec{v}_r \times \tilde{B}_\phi = \tilde{E}_z, \vec{E}_{MSE}$ is mostly in the \hat{z} direction $\rightarrow \widetilde{\Delta\lambda}_{Stark} \propto \tilde{B}_\phi$
 - $\vec{v}_r \times \tilde{B}_z = \tilde{E}_r$, changes angle of $\vec{E}_{MSE} \rightarrow$ measure polarization intensity
- Typically for tokamaks $\tilde{B}/B \sim 10^{-5}$, 100x smaller than broadband \tilde{E}/E
- However for EHO, $\tilde{B}/B \sim 10^{-4}$, coherent, low frequency (~ 10 kHz)
- \tilde{B}/B and \tilde{E}/E may be distinguishable using different beam energy components
 - $\Delta\lambda_1 \propto v_b \tilde{B} + \tilde{E}_{int}, \Delta\lambda_2 \propto \frac{1}{2} v_b \tilde{B} + \tilde{E}_{int}$




Diagnostic Requirements



\tilde{E}/\vec{E}_{MSE} in Fusion Grade Plasmas is $\sim 10^{-3}$

Experiment	$T_{e,0}$ (keV)	B (T)	a (m)	\tilde{E}/E_{MSE}
NSTX-U	$\sim 2-4$ (?)	1	0.6	$1 - 2 \times 10^{-3}$
DIII-D	2-5	2	0.7	$0.5 - 1 \times 10^{-3}$
Pegasus	~ 0.3	0.3	0.35	$0.7 - 1 \times 10^{-3}$

- Tokamak drift wave turbulence scaling gives $\tilde{E} \propto T_e/a$
- \tilde{E} turbulence broadband, majority of fluctuation power < 300 kHz
- \tilde{E} , \tilde{n} rise from core to edge
- \tilde{E}/E_{MSE} at or below photon noise floor for BES (typical rms noise $\sim 0.1\%$)
 - Two independent but spatially correlated measurements (e.g. $\langle \tilde{E}\tilde{n} \rangle$) made simultaneously can suppress incoherent photon noise another $\sim 10\times$



C. P. Ritz, et. al., Phys. Rev. Lett., **62**, 1989.



$\Delta\tilde{\lambda}_\pi$ Spectrometer Requirements are Formidable

- Resolution:

- Need ~ 8 spectral bins to resolve 2 gaussian-like π components
- $\Delta\lambda_{H\alpha}^\pi \sim 6 \text{ \AA}$ in D3D giving approximate spectral resolution of 0.75 \AA , $R = \frac{\lambda}{\delta\lambda} \sim 9000$

- Throughput:

- Matched to collection optics, $U = 0.016 \text{ cm}^2\text{sr}$ per D3D BES spatial point
- 2 spatial points desired


- Compatible detector system:

- $\sim 500 \text{ kHz}$ time response

- Mitigation of sightline-DNB geometric broadening



- Spectrometer optics conjugate to plasma collection lens \rightarrow coherent fiber bundle

D3D BES Spectrum:



Spatial Heterodyne Spectroscopy Technique Utilized for $\Delta\tilde{\lambda}_\pi(t)$ Measurement


- Self scanned, 2 beam interferometer
- Input wavelengths heterodyned around Littrow wavelength
- Field widening prisms increase throughput >10x




$$I(x) = \int_0^\sigma B(\sigma)(1 + \cos[2\pi(4|\sigma - \sigma_0|x \tan \theta_L)])d\sigma$$

Spatial Heterodyne Spectroscopy Achieves High Resolution and Throughput

- Phase I SHS design points:
 - $R = \sigma/\delta\sigma \approx 2Wd, U = 2\pi\eta A/R$
 - For \tilde{E} : Need approximately $R \sim 9000$
 - Design studies of SHS indicate low groove density grating maximizes sensitivity to \tilde{E}/E_{MSE}
 - Phase I parameters: grating width $W = 75$ mm and 50 g/mm, $R \sim 7500$ @ ~ 0.05 cm²sr
- Future design to utilize field widening prisms
 - Increases $U \sim 100\times$ at the same resolving power
 - Smaller SHS and multiple spatial points



Compact monolithic SHS design scales favorably to multiple plasma spatial points:



SHIMMER SHS: $R \approx 2.5 \times 10^4$, $U \approx 0.1$ cm²sr

Harlander, J. M., Roesler, F. L., Englert, C. R., et al. (2003).
Applied Optics, 42(15), 2829–2834.

Phase I Progress




Phase I Optical Layout

At tokamak:


Fiber run to D3D BES room:

Spectrometer and detector:



- Planned full system deployment to D3D for late 2017
- Goal is to validate spectrometer design with measurement of low frequency (~ 10 kHz) electric or magnetic field turbulence


Phase I deployed to D3D for Evaluation in Tokamak Environment





Initial Performance Characterization

- Single input wavelength leads to single spatial frequency
- Field lens corrects intensity vignetting and fringe distortion due to demagnification lens
- Hydrogen-Deuterium Giessler tube light source to be used for future resolution validation
 - Separation of H_{α} (656.279 nm) and D_{α} (656.1 nm) ~ 0.18 nm is close to the desired resolution for the measurement

Without field lens





With field lens



Broadband spectral input leads to shrinking of interferogram


- Variable line width put into SHS using tunable light source
- Spectrally broad input effects interferogram envelope
 - $I(x) = e^{-wx^2} [1 + \cos(8\pi(\sigma - \sigma_L))x \tan \theta_L]$
 - Multiple broad lines has addition effect
- For fluctuation measurement, desire every pixel in detector to provide meaningful information
 - Design studies of spectrometer sensitivity push groove density down




Compensation Plate Needed to Maintain Fringe Contrast, Match Fiber Etendue

- Spectrometer designed to match single fiber etendue
 - $U_{BES}=0.016 \text{ cm}^2\text{sr}$
 - Full etendue SHS aperture size $\sim 4.5 \text{ mm}$ dia. imaged by 75 mm dia. $f/2$ collimator
- Compensation plate required to maintain fringe contrast due to large aperture
- Plate installed in front of 50:50 beam splitter
 - Thickness equal to beam splitter plate, anti-reflection coated, $\lambda/4$ flatness


Small aperture



Full aperture



Fringe contrast with compensation plate:



D3D BES Lab vibration environment

- D3D BES lab elevated, thin walls, near HVAC equipment
- Active vibration isolation table nullifies table resonance and vibrations less than 20 Hz
- Further damping of horizontal vibrations may be necessary

SHS central fringe movement due to vibrations




Table sensors: vertical








Table sensors: horizontal



Sharp cutoff D_α filter required to isolate full energy components



- Filter designed by Alluxa
 - 90-10% cutoff < 0.25 nm
 - $> 90\%$ transmission in passband
 - 3 in diameter
- Allows for selection of only full energy beam components (improves sensitivity) or all components (distinguish between \tilde{B} and \tilde{E})
- Movement of passband over large filter area shallows cutoff
 - 0.4 nm 90-10% cutoff useable, misses design point
 - Working with manufacturer on fix



First Light Into Phase I SHS

- First light into SHS using fiber channel next to BES fiber array
 - D3D run for two days in July
- Observed 150L beam modulation
- SHS noise level estimated to be ~0.5%
 - $\sim 4 \times 10^7$ photo- e^- in 2 ms $\rightarrow 2 \times 10^4$ photo- $e^-/\mu\text{s}$
 - At 300 kHz BW, $\frac{d\tilde{N}}{dt} = \sqrt{2 \frac{dN}{dt} BW}$ gives $d\tilde{N}/dt = 346$ or ~0.5% rms noise
- RMS noise level likely to improve with filter and compensation plate fix


SHS central fringe intensity:



Geometric Broadening Compensation




Geometric Broadening Limits Diagnostic Sensitivity to \tilde{E}




- λ shift due to beam-collection optic viewing geometry is linear across window
 - Shift of ~ 0.3 nm (~ 7.5 cm $^{-1}$) across window
- Can be removed numerically or physically in spectrometer
- Both removal techniques require collection optic be imaged to a place inside spectrometer


Geometric broadening on DIII-D:



Sensitivity to \tilde{E}/E_{MSE} at different window sizes:



Collection Window Conjugate to Input Aperture Required to Remove Geometric Broadening



- 5 mm aperture size sets plasma collection volume (~2 cm)
- Fore optics image the plasma to the 5 mm aperture while imaging the BES collection lens to a ~10x10 coherent fiber bundle of approximately the same size as the current BES fiber bundles (~4.5 mm diameter) → reduce geometric broadening by ~10x


Geometric Broadening Compensation Achieved by Offsetting SHS Aperture

- Shift in real wavelength across BES collection lens can be negated by offsetting aperture (standard practice in atmospheric wind Fabry-Perot spectrometers see Trauger *et al.*)
- At spectrometer, window image rotated 90° so that $\Delta\sigma_{geo}$ is perpendicular to dispersion plane ($\Delta\sigma_{geo}(\alpha)$, where α is angle off dispersion plane)
- Grating equation off-axis has groove spacing (a) change proportional to α ,

$$\sigma a \cos \alpha (\sin \theta_{in} + \sin \theta_{out}) = m$$
- Solving for wavenumber shift while keeping OPD phase constant:


$$\delta(\Delta\sigma) = \frac{\Delta\sigma_0}{\cos^2(\theta_L)} \tan \alpha_0 \delta(\alpha)$$

- $\delta(\alpha)$ becomes the field of view of the spectrometer β



Doppler compensation concept in SHS:

~5x5 coherent fiber bundle



J. T. Trauger., Appl. Opt., **11**, 1972.

Summary: Moving Towards Diagnostic For \tilde{E} and \tilde{B} in High Temperature Magnetically Confined Plasmas Using Spatial Heterodyne Spectroscopy

- \tilde{E} measurements appear feasible using high speed measurements of motional stark multiplet
 - Horizontal view gives \tilde{E}_z , off midplane gives \tilde{E}_R
 - \tilde{B} measurement also possible using Stark multiplet
- Phase I SHS design completed, near deployment for first fluctuation data
 - First light through phase I SHS allowed testing of integrated system in tokamak environment
 - Need for improved filter transmission, compensation plate → changes in progress
- Next steps:
 - Validate technique with low frequency (~ 10 kHz) \tilde{E} or \tilde{B} measurement at D3D
 - Design Phase 2 spectrometer with geometric broadening compensation, run single coherent fiber bundle to D3D BES diagnostic lab



Reprints

*Reprints of this and other PEGASUS presentations
are available online at*

http://pegasus.ep.wisc.edu/Technical_Reports

