ABSTRACT



A CCD imaging system is currently being developed for $T_e(0,t)$ and bolometric measurements on the Pegasus Toroidal Experiment. Soft X-rays (E<10 keV) can be directly measured using a back-illuminated, thinned CCD to count individual photons and determine their energies via pulse-height analysis. Time resolution is obtained by moving exposed pixels across the CCD array behind a masked-off area on the chip. The first tests are done with a 512x512 array giving four spatial points and dt=2.5 msec. The mask is designed to optimize the tradeoff between time, energy, and spatial resolution. The presently available camera system is vacuum-coupled to the machine behind a 4 mil Be filter looking through a pinhole of diameter 0.16mm allowing measurements between R=15 and 45cm in the plasma. This same design can be extended to give a 1-D bolometric measurement by enlarging the pinhole and removing the Be filter. Initial measurements for both bolometry and electron temperature are in progress.

Outline



Concepts

- CCD-based SXR PHA on CHS/LHD
- Added capability of time evolution
- Bolometry application

SXR Modeling and T_e fits

- X-ray spectrum model and measurement simulations
- Noise Analysis
- Simulated data results

Hardware Implemetation

- Camera schematic and vacuum interface
- Lines of sight and machine view
- Initial tests of SXR sensitive CCD


Electron Temperature and Radiated Power Will Be Measured With a CCD Camera

- Soft x-ray sensitive CCD measures individual photons in pulse-height analysis mode
 - Provides SXR spectrum in the ~2-10 keV energy range
- The CCD/PHA is accurate and cost effective
 - Assuming a thermal electron distribution -fair assumption for PEGASUS
 - CCD/PHA validated in CHS/LHD work
- Bolometry measurements also available
 - Measure total energy flux with no spectral resolution to give time evolving radiated power measurement
- Provide critical measurement
 - $T_e(R,t)$ needed for stored energy, confinement
 - P_{rad} monitors impurity influence on plasma behavior

CCD/PHA Measurement of T_e Demonstrated on CHS/LHD



- CCD camera system implemented by Liang¹ et al. to measure T_e
 - Back-illuminated 1024 x 1024 CCD array directly detects individual soft x-rays
 - Beryllium filters and pinholes reduced the flux of photons to allow single photon detection
 - Multiple shots averaged to get adequate statistics
 - T_e obtained by fit to SXR continuum




Liang, Y., et al, Review of Scientific Instruments, vol. 71 no. 10 (Oct. 2000)

Pegasus Will Implement A CCD/PHA System With Time Resolution Capability



Implement partial masking to give time-resolved spectra with controlled frame transfer

- Use fast parallel line shift to move exposed areas behind mask
 - achieve ~1-2 ms time resolution



Requirements


- Sufficiently large number of pixels to provide decent photon statistics while avoiding pulse pile-up
- Maximize storage area to give large number of time intervals and spatial points
- Very low noise readout for good energy resolution

Radiated Power Measurements Made by Direct CCD Exposure Bolometer



 Direct exposure (i.e. no Be filter) of thinned CCD array has sensitivity equivalent to SXR diodes

 Sensitive to bulk of the impurity radiation



- Not restricted to avoidance of pulse pile-up
 - multiple time and spatial points


Full spatial profile of P_{rad} with compact camera


X-Ray Code Used to Model Pegasus Plasmas



 Initial SXR spectral simulations models the emissivity of Pegasus plasmas

Input model N_e(R), T_e(R), N_z/N_e






Calculate line integrated intensity as a function of energy and tangency major radius


Fit model spectrum to obtain Te

$$I \sim \sim \frac{1}{E\sqrt{T_e}}e^{-\frac{E}{\sqrt{T_e}}}$$

Deduction of SXR Spectrum Depends On Several Factors







- Photon flux reduced to appropriate level by apertures and filters
 - Be filter thickness adjusted to avoid impurity lines (<2 keV)
 - Pinhole area changed to reduce light flux and prevent pulse pileup
- Total number of pixels required for detection of N photons set by need to minimize pulse pileup
 - $N_{pixel} > N^2$ (Liang, et al.)

Slow-Scan CCD/PHA Decector Provides A Good Signal To Noise Ratio For Expected Pegasus Parameters

Noise contributions


- Dark current and noise will build up throughout the exposure to give noise ~10 electrons
- Photon noise is $\sim \overline{N}$ where N is the total number of photons collected
- Readout noise will be ~10 electrons

Total Noise


- Total noise (dark + readout) is ~15 electrons
- Photons are converted to electrons at ~3.2 eV/electron
- This results in a optimal spectral resolution of $E \sim 50 \text{ eV}$

Data Handling and Error Estimates Obtained With Simulations From X-Ray Code







A) Raw intensity at R_{tan}=30 cm



C) Random photon events distributed using part B as a probability distribution; N $<\sqrt{N_{pixels}}$; here N =180




B) Part A x QE x Filter transmission




D) Simulated spectrum with QE and transmission correction; fit to get T_e


Fits Yield Good Representation of Te






- •Fitting intensity to the local emissivity model gives a good representation of T_e at R_{tan} for $R > R_o$
- •Monte Carlo analysis gives T_e/T_e of 5-20%
 - Model determines number of photons required for a given error at a certain temperature





An Initial CCD/PHA System Has Been Assembled






- Mounted behind gate valve on Pegasus to allow access to filter and aperture
- First generation: Princeton Instruments 512 x 512 back-thinned CCD

Several Different Fields of View Are Available Depending On the Application






-fewer time points


-shot averaging

Mask Design Is A Key Element In Getting Time Resolution



- The CCD mask was made with two slots that expose two 32 x 512 pixel areas which are added together to get more pixels for improved photon statistics
- These exposed rows are shifted down during the plasma discharge to obtain a time record of data
- The parallel shift time is ~80 μs, which results in eight time points of t = 2.5 ms through the discharge
- Binned in the horizontal direction for a single spatial point
 - $N_{pixels} = 32768$ N max = 180 photons
- Pinhole diameters and locations are adjusted depending on the required photon flux and field of view





Summary



- Te and radiated power measurements will improve understanding of Pegasus plasmas
 - Time-resolved absolute measurement of T_e with CCD/PHA
 - $P_{rad}(R,t)$ to monitor impurity influence

- CCD/PHA has proven accurate
 - Current system will give $T_e(0,t)$ within ~10%
 - t~2.5 ms for eight time points through the discharge
 - T_e(R,t) (for R>R_o) shows good agreement with TS on CHS/LHD

- Testing of prototype system in progress
 - Shows sensitivity to x-rays from calibration source